

Children's Mercy Kansas City

SHARE @ Children's Mercy

Research Days

GME Research Days 2020

May 13th, 11:30 AM - 1:30 PM

Performance Analysis of Three Bioinformatic Variant Callers Using a Somatic Reference Standard

Binu Porath

Binu Porath

Let us know how access to this publication benefits you

Follow this and additional works at: <https://scholarlyexchange.childrensmercy.org/researchdays>



Part of the [Bioinformatics Commons](#), and the [Genomics Commons](#)

Porath, Binu and Porath, Binu, "Performance Analysis of Three Bioinformatic Variant Callers Using a Somatic Reference Standard" (2020). *Research Days*. 6.

https://scholarlyexchange.childrensmercy.org/researchdays/GME_Research_Days_2020/researchday3/6

This Poster Presentation is brought to you for free and open access by the Conferences and Events at SHARE @ Children's Mercy. It has been accepted for inclusion in Research Days by an authorized administrator of SHARE @ Children's Mercy. For more information, please contact hlsteel@cmh.edu.

Performance Analysis of Three Bioinformatic Variant Callers Using a Somatic Reference Standard

Binu Porath¹, Byunggil Yoo², Neil Miller², Midhat Farooqi¹

¹Department of Pathology and Laboratory Medicine, ²Center for Pediatric Genomic Medicine
 Children's Mercy Hospital, Kansas City, MO

Introduction

- Clinical next-generation sequencing assays require rigorous testing of bioinformatics pipelines to ensure that they faithfully detect relevant variants and minimize false positive and negatives.
- Here we compare the performance of 3 different bioinformatic variant callers for somatic variant detection (SVD) using a cell line-derived reference standard (RS).

Methods

Reference Standard (RS) characteristics

- OncoSpan Reference Standard (RS)
- Cell line-derived
- 386 variants in 152 genes
- Variants between 1-100% variant allele frequency (VAF)

Sequencing and alignment

- Illumina NovaSeq 6000
- Whole Genome Sequencing (WGS), 100X coverage
- Reads aligned to GRCh37/hg19 using Illumina DRAGEN Bio-IT (DB) platform

Analysis parameters

- Mostly default parameters
- Tumor-only analysis
- Turned off "contamination filter" to account for admixed RS

Variant calling on BAM file

- DRAGEN in Germline mode (DBG)
- DRAGEN in Somatic mode (DBS)
- GATK4 Mutect2 (MT2)

Performance Comparison




Figure 2

- DBG detected only 34/70 variants with a variant allele frequency (VAF) of more than 30%.
- DBG was removed from further analysis and we compared and contrasted DBS and MT2.

Filters Used

Panel of Normals (PoN)

- 82 normal samples used to account:
 - Misread bases
 - Sequencing artifacts
 - Misaligned reads
 - Position based filtering

Clustered events

- Filtered out variants if another variant was "nearby"- 15-100bp away

Bad haplotype/germline presence

- Filtered variants commonly present in healthy control databases
- gnomAD

Improper Filtering

DBS Improper Filtering

17 variants
 Clustered events

MT2 Improper Filtering

9 variants
 Germline risk





Figure 4

- MT2 inappropriately filtered out 9 variants because of "germline risk" even when the variant was absent in gnomAD.
- DBS improperly filtered out 17 variants due to "clustered events" if another variant was 'nearby' (anywhere from 15 to 100 base pairs away).
- Overall, they were concordant for 69% of the variant calls.

Summary & Conclusions



Results

OncoSpan Reference Standard (RS) variants




Figure 1: To focus on clinically relevant RS variants, we removed synonymous variants (n=268) and then any variant with a minor allele frequency (MAF) >1% in gnomAD (n=48). The remaining 70 variants were assessed.

Undetected/Filtered Variants

Breakdown of Undetected/Filtered Variants by DBS and MT2




Figure 3: Both DBS and MT2 did not detect 2 variants present at a VAF <3%, 6 variants 'present' in the PoN (filtered based on position, not nucleotide change), and appropriately filtered out 2 variants due to a bad haplotype and germline presence.

References

- Miller NA et al. A 26-hour system of highly sensitive whole genome sequencing for emergency management of genetic diseases. *Genome Med.* 2015; 7:100
- https://software.broadinstitute.org/gatk/documentation/tooldocs/4.beta.4/org_broadinstitute_hellbender_tools_walkers_mutect_Mutect2.php