Children's Mercy Kansas City

SHARE @ Children's Mercy

Manuscripts, Articles, Book Chapters and Other Papers

7-2018

Evaluation of the illumigene mycoplasma direct DNA amplification assay

Neena Kanwar Children's Mercy Hospital

Morgan A. Pence Cook Children's Medical Center

Donna Mayne Sacred Heart Hospital - Pensacola

Jeff Michael Children's Mercy Hospital

Rangaraj Selvarangan Children's Mercy Hospital

Let us know how access to this publication benefits you

Follow this and additional works at: https://scholarlyexchange.childrensmercy.org/papers

Part of the Medical Genetics Commons, and the Respiratory Tract Diseases Commons

Recommended Citation

Kanwar N, Pence MA, Mayne D, Michael J, Selvarangan R. Evaluation of the illumigene Mycoplasma Direct DNA Amplification Assay. J Clin Microbiol. 2018;56(7):e01930-17. Published 2018 Jun 25. doi:10.1128/JCM.01930-17

This Article is brought to you for free and open access by SHARE @ Children's Mercy. It has been accepted for inclusion in Manuscripts, Articles, Book Chapters and Other Papers by an authorized administrator of SHARE @ Children's Mercy. For more information, please contact hlsteel@cmh.edu.

BACTERIOLOGY



Evaluation of the *illumigene* Mycoplasma Direct DNA Amplification Assay

Neena Kanwar,^a ⁽ⁱ⁾ Morgan A. Pence,^b Donna Mayne,^c Jeffrey Michael,^a Rangaraj Selvarangan^a

^aChildren's Mercy Hospital and Clinics, Kansas City, Missouri, USA ^bCook Children's Medical Center, Fort Worth, Texas, USA ^cSacred Heart Hospital, Pensacola, Florida, USA

Journal of

SOCIETY FOR MICROBIOLOGY Clinical Microbiology®

AMERICAN SOCIETY FOR

ABSTRACT Mycoplasma pneumoniae is a common cause of community-acquired pneumonia. The illumigene Mycoplasma Direct (iMD) DNA amplification assay is a qualitative in vitro test utilizing loop-mediated isothermal amplification (LAMP) technology for the direct detection of *M. pneumoniae* DNA in respiratory specimens. The iMD assay does not require the preextraction of nucleic acids from specimens, which is a prerequisite step for the previously approved illumigene Mycoplasma (iM) assay. The aim of this prospective multicenter study was to evaluate the performance characteristics of the newly developed iMD assay, compared with the iM assay. Subjects with symptoms of upper respiratory illnesses suggesting M. pneumoniae infection were enrolled at three sites in the United States. Respiratory specimens were obtained using dual throat swabs. One swab was tested with the iMD assay at each enrollment site. Reference testing with the iM assay was performed by the manufacturer. Among 456 specimens tested, the iM reference method detected M. pneumoniae in 25 specimens (5.5%), while the iMD assay identified 34 specimens (7.5%) as M. pneumoniae positive. There were 10 false-positive results and 1 false-negative result with the iMD assay. The overall positive and negative agreement rates were 96.0% (95% confidence interval [CI], 80.5 to 99.3%) and 97.7% (95% CI, 95.8 to 98.7%), respectively. The overall agreement rate was determined to be 97.6% (95% Cl, 95.7 to 98.6%). We conclude that the iMD test results were comparable to the iM assay results. The removal of the DNA extraction step for the iMD assay simplifies testing, saves time, and reduces the costs of detecting M. pneumoniae from throat swabs, compared to the iM assay.

KEYWORDS Mycoplasma pneumoniae, illumigene assay, LAMP technology

Mycoplasma pneumoniae is a common cause of respiratory tract infections, especially in young adults and school-age children. Diagnosis of *M. pneumoniae* infection based on a patient's clinical presentation alone is difficult and unreliable (1). Prompt diagnosis is essential to initiate appropriate antibiotic therapy and infection control.

Several methods are utilized for *M. pneumoniae* detection. Culture-based methods are highly specific and sensitive but are time-consuming, and their sensitivity may vary depending on laboratory skills (2, 3). Serology-based assays (complement fixation assay, enzyme-linked immunoassay, enzyme-linked immunosorbent assay, and microparticle agglutination assay) are available. The limitations of serology-based assays include a lack of sensitivity due to a delayed antibody response to *M. pneumoniae* infection, a lack of response in older patients, highly prevalent background antibodies in healthy individuals, and cross-reactivity with other *Mycoplasma* species (2, 4, 5). Rapid antigen kits targeting specific proteins are also being utilized to diagnose *M. pneumoniae*. The sensitivity of rapid antigen assays has been reported to vary from 60% to 90% (2, 6, 7).

Received 12 December 2017 Returned for modification 10 January 2018 Accepted 19 April 2018

Accepted manuscript posted online 2 May 2018

Citation Kanwar N, Pence MA, Mayne D, Michael J, Selvarangan R. 2018. Evaluation of the *illumigene* Mycoplasma Direct DNA amplification assay. J Clin Microbiol 56:e01930-17. https://doi.org/10.1128/JCM .01930-17.

Editor Karen C. Carroll, Johns Hopkins University School of Medicine

Copyright © 2018 Kanwar et al. This is an open-access article distributed under the terms of the Creative Commons Attribution 4.0 International license.

Address correspondence to Rangaraj Selvarangan, rselvarangan@cmh.edu.

	No. of subj	M. pneumoniae		
Site no. and location	Total	Total positive	prevalence (%)	
Site 1, Florida	152	4	2.6	
Site 2, Texas	49	4	8.1	
Site 3, Missouri	255	26	10.2	
Total	456	34	7.5	

TABLE 1 Total subjects enrolled and *M. pneumoniae* prevalence across the three sites

Nucleic acid amplification testing (NAAT) assays have gained popularity due to increased assay sensitivity (1, 2), which enables detection at an early stage of infection. The loop-mediated isothermal amplification (LAMP) assay is a NAAT assay that is simple, easy to use (no thermocycler is needed), and highly sensitive, with a quick turnaround time (2). The illumigene Mycoplasma (iM) assay (Meridian Bioscience, Inc., Cincinnati, OH) utilizes LAMP technology for the detection of *M. pneumoniae* from throat and nasopharyngeal swab specimens. The primers of this diagnostic kit target a 208-bp DNA sequence found in the intracellular protease-like protein gene of the M. pneumoniae genome (1). A byproduct of the continuous isothermal amplification reaction is magnesium pyrophosphate. This white precipitate leads to turbidity in the reaction mixture. The illumipro-10 incubator/reader monitors the changes in the absorbance characteristics. Significant changes in the absorbance of the reaction mixture indicate the presence of the target gene. A study found the sensitivity and specificity of the iM assay to be 100% and 99%, respectively, compared with the culture method (1). Another recent study found both sensitivity and specificity of the iM assay to be 100%, compared with the FilmArray respiratory panel (bioMérieux, France) (8).

The company manufacturing the iM assay recently developed an improved version, the *illumigene* Mycoplasma Direct (iMD) assay, which does not require a Qiagen DNA extraction step. The gene target and assay chemistry are identical for the iM and iMD assays. This multicenter clinical trial compared the clinical performance of the iMD assay with that of the iM assay, which requires nucleic acid extraction. This is the first study evaluating the performance of the iMD assay.

MATERIALS AND METHODS

Study design. This prospective multicenter clinical trial was conducted at three different sites (Florida, Texas, and Missouri) across the United States. Patients with upper respiratory tract infections that might be attributed to *M. pneumoniae* were enrolled after informed consent was obtained. Fresh throat swab specimens from both male and female patients were collected between August 2015 and January 2016. This study was designed to evaluate the performance of the newly improved iMD assay, compared with the iM assay. The study protocols were approved by each hospital's institutional review board.

Specimens. Specimen inclusion criteria included (i) specimens collected from patients with symptoms of upper respiratory tract infections that might be attributed to *M. pneumoniae* or patients suspected of having *M. pneumoniae* infection, (ii) two throat swabs from each subject enrolled, one transported in nonnutrient transport medium and other in M4 medium, and (iii) written informed consent provided by the subject. Specimen exclusion criteria included (i) specimens received in the laboratory in unsatisfactory containers or conditions, (ii) multiple sets of specimens from the same patient at different office visits, (iii) specimens received in the laboratory with less than two throat swabs per subject, and (iv) specimens stored in unapproved transport medium types or combinations. Demographic data collected at enrollment included age, gender, date of symptom onset, antibiotic use in past 4 weeks, and other medication use. No other clinical or radiological findings were recorded.

A total of 456 subjects at the three different sites were included in the study. A summary of the subjects enrolled and *M. pneumoniae* prevalence across the three sites is provided in Table 1. Patients 3 weeks to 97 years of age with upper respiratory tract illnesses that might be attributable to *M. pneumoniae* were included. A dual-swab (rayon tip) collection device was used to collect one specimen per patient. One swab was stored in liquid Amies medium (catalogue no. 220105; Becton, Dickinson and Company, Franklin Lakes, NJ, USA) and the second was stored in M4 transport medium (catalogue no. 12500; Thermo Fisher Scientific; Remel Products, Lenexa, KS, USA). Specimens were stored at 2 to 8°C prior to testing. At each study site, specimens were deidentified and assigned a unique, site-specific, study identifier. Testing was performed with the iMD assay with one throat swab at each study site. Operators at each study site tested external positive- and negative-control samples every day prior to performing the iMD assay with the study samples. The second swab, in M4 medium, was shipped to the

TABLE 2 *M. pneumoniae* prevalence according to age and gender across all three sites, as determined with the iMD assay

	No. of sub					
Age and gender	Site 1 (Florida)	Site 2 (Texas)	Site 3 (Missouri)	Total	Total positive	<i>M. pneumoniae</i> prevalence (%)
Age group						
0–1 mo	0	0	7	7	0	0.0
2 mo to 2 yr	9	12	136	157	5	3.2
3–12 yr	20	28	99	147	22	15.0
13–21 yr	16	9	13	38	5	13.2
22–65 yr	86	0	0	86	2	2.3
>65 yr	21	0	0	21	0	0.0
Gender						
Male	63	24	139	226	18	8.0
Female	89	25	116	230	16	7.0
Total	152	49	255	456	34	7.5

manufacturer for testing with the iM assay. The results from the two assays were compared to determine the performance of the iMD assay.

illumigene Mycoplasma Direct assay. The swab stored in red-capped nonnutrient liquid Amies transport medium was utilized to perform the iMD assay, at the clinical site, within 72 h after sample collection. Briefly, the swab was inserted by breaking the swab handle in the SMP Prep tube provided with the kit. The SMP Prep tube was vortex-mixed for 10 s. Five to 10 drops from the SMP Prep tube were squeezed into a 1.5-ml heat treatment tube. The heat treatment tube was heated at 95°C for 10 min, followed by 10 s of vortex-mixing. Fifty microliters of this heat-treated sample was transferred to the test chamber as well as the control chamber of the *illumigene* test device. The device was closed, and all air bubbles were removed by gently tapping the test device. The *illumigene* test device was inserted into an *illumipro*-10 incubator/reader for direct detection of *M. pneumoniae*. Results were displayed on the instrument as positive or negative at the conclusion of the run, in less than 1 hour.

illumigene Mycoplasma assay. Specimens stored in M4 medium were shipped (transport temperature, 2 to 8°C) to the company within 4 days after sample collection and were tested within 14 days. The Qiagen QIAamp DSP DNA minikit was utilized to extract DNA from 150 μ l of the specimen. The comparator iM assay was performed by the company according to the assay package insert.

Statistical analysis. (i) Descriptive statistics. Overall characteristics of the subjects enrolled and the *M. pneumoniae* prevalence at each individual site were determined with reference to the iMD assay. *M. pneumoniae* prevalence was also determined according to age and gender for each of the three sites.

(ii) Analytical statistics. Data tables (two by two) were utilized to determine the rates of positive and negative agreement between the iM assay and the iMD assay. Analysis of the performance characteristics (positive and negative agreement rates), with 95% confidence intervals (CIs), was performed using the Vassarstats website (http://vassarstats.net/clin1.html).

RESULTS

A total of 456 subjects with symptoms of upper respiratory tract illness, from three different sites, were enrolled in the study. The highest prevalence was observed at the Missouri site (10.2%), followed by the Texas site (8.1%) (Table 1). There was an approximately equal distribution of male and female subjects enrolled in the study (female, n = 230 [50.4%]; male, n = 226 [49.6%]) (Table 2). The median age of the patients included in the study was 5 years (range, 0.06 to 97 years). The highest *M. pneumoniae* prevalence was observed in the age group of 3 to 12 years (22/147 patients [15%]), compared to the rest of the study population (P < 0.001). The association of age with *M. pneumoniae* prevalence has been documented in previous studies (9–11). *M. pneumoniae* prevalence rates among female and male patients (7.0% and 8.0%, respectively) were found to be comparable (P = 0.7). Detailed age- and gender-specific prevalence rates across all sites are listed in Table 2.

Among the 456 specimens tested, the iMD and iM assays detected *M. pneumoniae* in 34 specimens (7.5%) and 25 specimens (5.5%), respectively. There were 10 false-positive specimens and 1 false-negative specimen detected by the iMD assay. The rates of positive and negative agreement between the iMD assay and the iM assay were 96.0% and 97.7%, respectively. Discrepancy analysis was performed by repeat testing of a second aliquot of the discrepant sample with the iM assay. Four of 10 false-positive

iMD assay result	No. of samples ^a					
	Positive iM assay result	Negative iM assay result	Total			
Positive	24	10 ^b	34			
Negative	1 ^c	421	422			
Total	25	431	456			

TABLE 3 Performance characteristics of the iMD assay, compared with the iM assay (with Qiagen extraction)

^oOverall agreement, 97.6% (95% Cl, 95.7 to 98.6%); positive agreement, 96.0% (80.5 to 99.3%); negative agreement, 97.7% (95.8 to 98.7%).

^bFour of 10 samples were identified as positive by the iM assay after testing with an additional frozen swab collected from the patients for discrepancy analysis.

cRepeat testing with the iM assay, with the original patient sample and an additional frozen sample, produced negative results.

specimens were identified as positive by the iM assay after retesting with an additional frozen sample. Repeat testing of 1 false-negative specimen with the iM assay, using the original patient sample and an additional frozen sample, produced negative results in both cases.

Detailed performance parameters for all study subjects and specific parameters for the study population are presented in Table 3. No invalid runs were observed during the study. Results of standard-of-care tests for *M. pneumoniae* detection (an enzyme immunoassay and PCR for the first site, the iM assay for the second site, and the BioFire FilmArray respiratory panel for the third site) were available for 70 of the 456 subjects enrolled in the study, with 11 positive detections.

DISCUSSION

The objective of this multisite clinical trial was to evaluate the performance characteristics of the iMD assay, compared to the iM assay. This is the first study evaluating the performance of the iMD assay.

A previous study (8) in the United States compared the iM assay, the Prodesse ProPneumo-1 assay (Hologic Gen-Probe, San Diego, CA, USA), and the *Mycoplasma pneumoniae* P1 LightMix kit (TIB Molbiol, Howell, NJ, USA) with the FilmArray respiratory panel (bioMérieux, France) for detection of *M. pneumoniae* from pediatric clinical nasopharyngeal swab specimens. All of the three commercially available NAAT assays had similar sensitivities, specificities (100%, 100%, and 96% for the iM assay, the Prodesse assay, and the LightMix kit, respectively), and hands-on times. Another study compared the iM assay with culture, using frozen respiratory specimens from adults and children for whom historic culture results were available (1). The sensitivity and specificity of the iM assay after discrepancy analysis were observed to be 100% and 99%, respectively. In our study, the overall rate of agreement between the iMD assay and the iM assay was found to be 97.6%. Given the study results, it is expected that the iMD assay will perform comparably, compared with the existing diagnostic platforms.

Current FDA-cleared sample-to-answer molecular methods for *M. pneumoniae* detection are the nested multiplex PCR panel by BioFire Inc. and standalone LAMP assays (iM and iMD assays) by Meridian Bioscience, Inc. The iMD assay targets the intracellular protease-like protein gene by utilizing LAMP technology (no thermocycler is needed) and provides results in less than 1 hour. By comparison, the FilmArray respiratory panel targets the *tox* gene for *M. pneumoniae* detection. Advantages of the iMD assay over the FilmArray respiratory panel for *M. pneumoniae* detection are that both instrument costs and costs per test are lower for the iMD test and the iMD assay has a higher throughput (it can process 10 specimens, compared with a single specimen, at one time), making it ideal for testing in an outpatient setting. One study evaluated the performance of the FilmArray respiratory panel for detection of *M. pneumoniae* versus a laboratory-developed real-time TaqMan PCR assay targeting the *p1* gene (12) and determined the positive agreement rate to be 100% (95% CI, 70.1 to 100%), which was comparable to the value of 96.0% (95% CI, 80.5 to 99.3%) for both iM and iMD assays obtained in our study. We evaluated 25 true-positive *M. pneumoniae* isolates, compared with 9 positive isolates utilized in the FilmArray respiratory panel study (12); this resulted in relatively more robust estimates, as indicated by the 95% CIs.

Several studies conducted in Europe have compared various commercially available NAAT assays for *M. pneumoniae* detection (13–15), and a summary of commercially available NAAT assays for detection of *M. pneumoniae* was provided in a recent minireview article (16). These NAAT assays need specialized trained personnel, expensive setups with thermal cycling amplification platforms, and a molecular biology facility. In contrast, the LAMP technology utilized in the iMD assay allows specific and continuous DNA amplification under isothermal conditions. The iMD assay is a sensitive assay that does not require an additional DNA extraction step, making it faster and more economical than other NAAT assays that require the DNA extraction step. Overall, the iMD assay is a simple, convenient, and rapid molecular assay for detection of *M. pneumoniae* from throat swab specimens.

ACKNOWLEDGMENTS

Funding for this study was provided by Meridian Bioscience, Inc. We acknowledge the dedicated members of the enrollment team and laboratory technicians at each site for helping with this study.

REFERENCES

- Ratliff AE, Duffy LB, Waites KB. 2014. Comparison of the *illumi*gene Mycoplasma DNA amplification assay and culture for detection of *Mycoplasma pneumoniae*. J Clin Microbiol 52:1060–1063. https://doi.org/10 .1128/JCM.02913-13.
- Parrott GL, Kinjo T, Fujita J. 2016. A compendium for *Mycoplasma pneu-moniae*. Front Microbiol 7:513. https://doi.org/10.3389/fmicb.2016.00513.
- Diaz MH, Winchell JM. 2016. The evolution of advanced molecular diagnostics for the detection and characterization of *Mycoplasma pneumoniae*. Front Microbiol 7:232. https://doi.org/10.3389/fmicb.2016.00232.
- Youn Y-S, Lee K-Y. 2012. Mycoplasma pneumoniae pneumonia in children. Korean J Pediatr 55:42–47. https://doi.org/10.3345/kjp.2012.55 .2.42.
- Waites KB, Balish MF, Atkinson TP. 2008. New insights into the pathogenesis and detection of *Mycoplasma pneumoniae* infections. Future Microbiol 3:635–648. https://doi.org/10.2217/17460913.3.6.635.
- Miyashita N, Kawai Y, Tanaka T, Akaike H, Teranishi H, Wakabayashi T, Nakano T, Ouchi K, Okimoto N. 2015. Diagnostic sensitivity of a rapid antigen test for the detection of *Mycoplasma pneumoniae*: comparison with real-time PCR. J Infect Chemother 21:473–475. https://doi.org/10 .1016/j.jiac.2015.02.007.
- Yamazaki T, Kuroki H, Itagaki T, Iwata S, Tateda K. 2015. Evaluation of a rapid antigen detection kit targeting L7/L12 ribosomal protein for *Mycoplasma pneumoniae*. Kansenshogaku Zasshi 89:394–399. https://doi .org/10.11150/kansenshogakuzasshi.89.394.
- Chou RC, Zheng X. 2016. A comparison of molecular assays for *Mycoplasma pneumoniae* in pediatric patients. Diagn Microbiol Infect Dis 85:6–8. https://doi.org/10.1016/j.diagmicrobio.2015.12.013.
- 9. Chen ZR, Yan YD, Wang YQ, Zhu H, Shao XJ, Xu J, Ji W. 2013. Epidemi-

ology of community-acquired *Mycoplasma pneumoniae* respiratory tract infections among hospitalized Chinese children, including relationships with meteorological factors. Hippokratia 17:20–26.

- Othman N, Isaacs D, Kesson A. 2005. *Mycoplasma pneumoniae* infections in Australian children. J Paediatr Child Health 41:671–676. https://doi .org/10.1111/j.1440-1754.2005.00757.x.
- Defilippi A, Silvestri M, Tacchella A, Giacchino R, Melioli G, Di Marco E, Cirillo C, Di Pietro P, Rossi GA. 2008. Epidemiology and clinical features of *Mycoplasma pneumoniae* infection in children. Respir Med 102: 1762–1768. https://doi.org/10.1016/j.rmed.2008.06.022.
- Pierce VM, Elkan M, Leet M, McGowan KL, Hodinka RL. 2012. Comparison of the Idaho Technology FilmArray system to real-time PCR for detection of respiratory pathogens in children. J Clin Microbiol 50:364–371. https://doi.org/10.1128/JCM.05996-11.
- Touati A, Benard A, Hassen AB, Bébéar CM, Pereyre S. 2009. Evaluation of five commercial real-time PCR assays for detection of *Mycoplasma pneumoniae* in respiratory tract specimens. J Clin Microbiol 47:2269–2271. https://doi.org/10.1128/JCM.00326-09.
- Dumke R, Jacobs E. 2014. Evaluation of five real-time PCR assays for detection of *Mycoplasma pneumoniae*. J Clin Microbiol 52:4078–4081. https://doi.org/10.1128/JCM.02048-14.
- Pillet S, Lardeux M, Dina J, Grattard F, Verhoeven P, Le Goff J, Vabret A, Pozzetto B. 2013. Comparative evaluation of six commercialized multiplex PCR kits for the diagnosis of respiratory infections. PLoS One 8:e72174. https://doi.org/10.1371/journal.pone.0072174.
- Loens K, leven M. 2016. Mycoplasma pneumoniae: current knowledge on nucleic acid amplification techniques and serological diagnostics. Front Microbiol 7:448. https://doi.org/10.3389/fmicb.2016.00448.